

2018 Show SPECIAL

- Fall Prevention & Therapy Support Specialists
- Simple Solutions Designed to Keep Patients Safe
- Patented Heel Off™ Technology
- Low Cost and Clinically Effective Products

On all Fall Prevention & Therapy Support Capital Purchases

Therapy Support

OrthoCare Air Plus

Low Air Loss
 Alternating Pressure
 Heel Off™ Technology
\$599.00



Air Mattress Features

- Modular air cells, easy to change and clean.
- Medical grade TPU material, environmentally friendly, waterproof, & anti-bacterial
- Air ventilation hole design, blowing by tiny holes, moisture control, patient can feel more comfortable
- Double CPR straps quickly deflate the mattress for cardiac resuscitation first aid
- Integrated Safe and Secure Safety Cover, helps reduce the risk of falls.
- Whisper quiet pump can be programmed specifically for the patient's weight.
- System lock: helps prevent accidental adjustments

OrthoCare Air

Low Air Loss
 Alternating Pressure
\$599.00



Safe & Secure™ Safety Cover

\$100

- MACHINE WASHABLE
- FLUID-RESISTANT
- BREATHABLE- STRETCH MATERIAL REDUCES ODOR
- PREVENTS SHEARING
- SAFETY COVER
- 4 PADDED BOLSTERS
- 1 FOOT PILLOW
- 4 ATTACHED STRAPS WHICH ANCHORS TO ANY BED



Fall Prevention

\$90.00

Safe & Secure™ Safety Coverlet

All the features of the Safety Cover designed to operate with most air mattresses

Bolsters slide out for a Safe & Comfortable Egress



Questions or Comments?

CALL (888) 616-9811

Email: customerservice@orthotecmedical.com

FOR MORE INFORMATION VISIT: WWW.ORTHOTECMEDICAL.COM

Ask about a 14 Day No Risk
FREE Evaluation